

CASCADE 20L BEAD BATH

No More Contamination or Evaporation

WATERLESS LABORATORY BEAD BATH

IMPROVES TEMPERATURE
PERFORMANCE & SAVES ENERGY




LABS ARE BETTER WATERLESS

Cascade Sciences' 20L Bead Bath is a waterless laboratory bead bath that reduces contamination and the messy condensation associated with water baths. Bead baths don't evaporate, which improves temperature performance and saves energy.

FEATURES

- 15L Lab Armor Beads included
- High Thermal Efficiency
- Microbial/Biofilm Resistant
- Cross-Contamination Resistant
- Minimal Cleaning
- No Refilling
- No Rack or Stand Accessories Needed
- Angled Incubation Possible
- Constant Volume – No Evaporation
- Non-Corrosive
- Compatible with non-water tight vessels
- Independent build-in over-temp, burn-out protection
- Mobile or Field Use Compatible
- Non-Combustible
- NRTL Certified



UPGRADE YOUR WATER BATH WITH LAB ARMOUR BEADS



503-847-9047

info@cascadesciences.com
cascadesciences.com

VISIT OUR SHOWROOM!

6725 NE Evergreen Pkwy, Ste 106
Hillsboro, OR 97124

**UPGRADE YOUR
WATER BATH WITH
LAB ARMOUR BEADS**



2L THERMAL BEADS



4L THERMAL BEADS



8L THERMAL BEADS


Cascade 20L Bead Bath

SPECIFICATIONS


20L BEAD BATH

Temperature Range	Ambient +5°C to 80°C
Temperature Uniformity	±1.0°C at 37°C
Bead Capacity	15L Thermal Beads Included
Bath Dimensions (without lid)	
Height	8.5" (21.6cm)
Width	15.9" (40.4cm)
Depth	23.25" (59cm)
Tub Dimensions	6"W x 12"D x 17"H
Power Options	100V – 120V / 1ph / -6AMP 50/60Hz
	220V – 240V / 1ph / -3AMP 50/60Hz

Front View



Back View



503-847-9047

info@cascadesciences.com
cascadesciences.com

VISIT OUR SHOWROOM!

6725 NE Evergreen Pkwy, Ste 106
Hillsboro, OR 97124