

# LA SERIES BENCHTOP

BENCHTOP FUME HOOD

**PRECISION**  
EXTRACTION SOLUTIONS

By AGRIFY™

**16G SUPERSTRUCTURE**  
REMOVABLE SIDE PANELS



## LA SERIES BENCHTOP

BENCHTOP FUME HOOD

Precision's benchtop fume hoods provide the safest work environment possible with down to 60 ft/min consistent laminar airflow across the face of the hood. Available in three standard sizes and fully customizable in width, height, and depth to fit your application or space, these fume hoods will provide you with versatility, durability, and quality. All options are built with 16-gauge galvanized steel-pan construction, angled front post design, and removable side panels and come equipped with a chemical-resistant poly resin liner. All benchtop fume hoods are proudly made in the USA and proudly assembled in Michigan!

### SYSTEM HIGHLIGHTS



**PROUDLY MADE IN THE USA**  
Assembled in Michigan



**CHEMICAL RESISTANT FINISH**  
Powder coat & interior poly resin



**16G COLD-ROLLED STEEL**  
Fully customizable dimensions



**ASHRAE-110**  
Tested & approved

# LA SERIES BENCHTOP

## BENCHTOP FUME HOOD

### Features

- Fully customizable width, height & depth
- 16g cold-rolled steel superstructure
- Stainless steel sash frame & guides
- Angled posts for smooth airflow transition
- Three-piece baffle system with dampers
- LED explosion-proof lighting for improved visibility
- Chemical resistant poly resin liner
- Chemical resistant powder coat finish
- ASHRAE-110 tested and approved

### Optional Equipment

- Remote- controlled plumbing fixtures with chrome handles, pre-piped with flex hose or hard copper.
- Multiple sash options & stops
- Custom pass-through openings upon request
- Air/ Monitor Alarm
- Ceiling Enclosure
- Base Cabinet
- Blower

### Equipment Specifications

Model 4854	
<b>Exterior Dimensions</b>	48" L x 33.5" W x 54" H
<b>Worksurface</b>	48" L x 30" W x 1.25" H
<b>Hood Opening</b>	40" W x 30" T
<b>Hood Weight</b>	459 lb
<b>CFM (approx.)</b>	500*

Model 6054	
<b>Exterior Dimensions</b>	60" L x 33.5" W x 54" H
<b>Worksurface</b>	60" L x 30" W x 1.25" H
<b>Hood Opening</b>	52" W x 30" T
<b>Hood Weight</b>	497 lb
<b>CFM (approx.)</b>	650*

<b>Color Options</b>	Custom Color-Match Available
<b>Warranty**</b>	One-Year

Model 7254	
<b>Exterior Dimensions</b>	72" L x 33.5" W x 54" H
<b>Worksurface</b>	72" L x 30" W x 1.25" H
<b>Hood Opening</b>	64" W x 30" T
<b>Hood Weight</b>	547 lb
<b>CFM (approx.)</b>	800*

### Compliance

<ul style="list-style-type: none"> <li>• UL Listed 1805 - iQ Labs OGOY.MH28517</li> </ul>	<ul style="list-style-type: none"> <li>• ASHRAE-110 Tested &amp; Approved</li> </ul>
---	--

