

By AGRIFY

OIL SEALED ROTARY VANE VACUUM PUMP






LAB SOCIETY

ELITE RELIABILITY AND PERFORMANCE

OPTIMIZED FOR OPERATION



MODEL			LS-RV-12*	LS-RV-18*
Pumping Rate	m ³ /hr (cfm)	50 Hz	9.9 (5.8)	14.4 (8.5)
		60 Hz	12 (7.1)	17.4 (10.2)
Ultimate Pressure	Torr	Gas Ballast Closed	4 x 10 ⁻³	
		Gas Ballast Open	4 x 10 ⁻²	
Motor Power	kW (HP)		0.4 (0.5)	0.55 (0.75)
Voltage	V	Three Phase	N/A	
		Single Phase	110 ~220 V	
Oil Quantity	Liters		0.55 ~ 1.1	0.65 ~ 1.2
Inlet	KF		25	
Outlet	KF		25	
Ambient Temp.	°F		41° ~ 104°	
Weight	lbs. (kg)		56 (25.4)	74 (34)

*Explosion Proof models available

MODEL			LS-RV-36*	LS-RV-50*	LS-RV-70*	LS-RV-110*
Pumping Rate	m ³ /hr (cfm)	50 Hz	30 (17.7)	40 (23.6)	60 (35.3)	90 (53)
		60 Hz	36 (21.2)	48 (28.3)	72 (42.4)	108 (63.6
Ultimate Pressure	Torr	Gas Ballast Closed	4 x 10 ⁻³			
		Gas Ballast Open	4 x 10 ⁻²			
Motor Power	kW (HP)		1.1 (1.5)	1.5 (2.0)	2.2 (3.0)	3.7 (5.0)
Voltage	V	Three Phase	220V	/ 230 ~ 460V		
	V	Single Phase	110V	V N/A		
Oil Quantity	Liters		1.2 ~	1.2 ~ 2.8 2.5 ~ 4.2		
Inlet	KF		40			
Outlet	KF		40			
Ambient Temp.	°F		41° ~ 104°			
Weight	lbs. (kg)		142 (64)	147 (67)	197 (89)	229 (104)

*Explosion Proof models available

MODEL			LS-RV-200	LS-RV-300
Pumping Rate	m ³ /hr (cfm) 60 Hz		196 (115.4)	306 (180)
Ultimate Pressure	Torr	Gas Ballast Closed	4 x 10 ⁻³	
		Gas Ballast Open	4 x 10 ⁻²	
Motor Power	ŀ	kW (HP)	5.6 (7.5)	7.5 (10)
Voltage	V	Three Phase	230 ~ 460V	
	v	Single Phase	N/A	
Oil Quantity	Liters	Max.	25	28
	Liters	Min.	20	23
Inlet	ISO		63	
Outlet	KF		40	
Cooling Water Req.	L/hr.	L/hr. Water Temp. 68°F		248
Weight		lbs. (kg)	513 (233)	568 (258)



BOOSTER PUMPS

	MODEL		LS-RB-120-3	LS-RB-330-3	LS-RB-600-3
Pumping Rate	cfm	50 Hz	58	164	294
		60 Hz	70	194	353
Max Inlet Pressure	Torr	50 Hz	9.00	9.00	9.75
- Indx IIIIct i i coodiic		60 Hz	6.95	6.95	8.25
Max Differential Pressure	Torr	50 Hz	30.00	30.00	54.75
		60 Hz	24.75	24.75	45.00
Ultimate Pressure	Torr		3x10 ⁻³		
Motor Power kW	HP		0.4 (0.5)	0.75 (1.0)	2.2 (3)
Motor Voltage	3-Stage		230/460V		
Motor Amp Rating	230V (460V)		1.6 (0.8)	2.8 (1.4)	7.2 (3.6)
Oil Quantity	Liters		0.4	0.8	1.6
Cooling Water	L/min		N/A	2.0	2.0
	PSI		N/A	14.5	14.5
	Temp (F°)		41-86		
Inlet			VG50	VG80	VG80
Outlet			VF50	VF80	VF80
Ambient Temp			41-104		
Weight lbs (kg)			67 (30)	113 (51)	177 (80)